



CONSTANT TEMPERATURE OIL / WATER BATHS

Used for high temperature indirect sample heating inside the bath and also for external circulation. Supplied with top SS lid

Temp. Range : 5°C above amb. To 300°C

Inner Chamber Sizes (Die pressed)

A. 300mm L x 200mm W x 200mm H,

Working space : 200L x 200W x 200 H

Outer size : 400 x 300 x 500

Volume : 12 Lit.

Weight : 20 Kgs

B. 525mm L x 320mm W x 200mm H


Working Space : 400 L x 320 W x 200 H

Volume : 32 lit.

Outer size : 625 x 420 x 500

Weight : 36 Kg.

C. Custom sizes also available



STANDARD FEATURES

1. Inner of SS 304 construction
2. PID based Temperature control system with pt-100 sensor.
3. TEFC motor based Pump/Stirrer
4. Special SS immersion heaters for operation with lower oil levels.

OPTIONALS

- A. Profile controller
- B. Circulator / pump
- C. Customised Sample Holding arrangement
- D. Diffuser Plate
- E. All SS Construction
- F. Liquid Media

