



General Purpose Laboratory Water Baths



Inner Chambers made of S.S.304 with outer of elegantly powder coated C.R. Sheet with antislip rubber legs and best quality thermostat, load and mains indicators and long life 'U' type immersion heaters. Also available are optional features like water levels, digital temp.controllers, die pressed single piece inner chambers, stirrers, external circulation pumps.

These models are available in various sizes as below and can also be custom built as per requirement and application.

- * 14" x 10" x 4" 6 Holes with concentric rings
- * 16" x 12" x 4" 12 Holes with concentric rings
- * 8" Dia x 6" Ht Cylindrical with concentric rings
- * 10" Dia x 6" Ht Cylindrical with concentric rings
- * 12" Dia x 6" Ht Cylindrical with concentric rings
- * 250mm L x 150mm W x 175mm D for 2 Racks
- * 350mm L x 250mm W x 175mm D for 4 Racks
- * 450mm L x 250mm W x 175mm D for 6 Racks



D-428, 10th Main, Peenya II Stage, P.I.A,
BANGALORE-58. INDIA.
PH : 091-080-28360799. Fax :23349239
E@mail: tic_ci80@hotmail.com;
tic_ci@yahoo.co.in