



When Reliability matters, only

EQUITRON®

SEROLOGICAL WATER BATH DIGITAL

Applications:

- ❖ Ideal for sample pre-treatment in QC, for Serological tests in Pathology and all thermostating tests.

EQUITRON Advantages:

- ❖ Fully stainless steel, including bottom, makes it GMP standard & total rust proof.
- ❖ Easy to set & use - direct digital setting & readout.

Construction:

- ❖ Double walled full SS construction with glasswool insulation.

Standard Features:

- ❖ Temperature Range: Ambient +5°C to 90°C.
- ❖ Digital Temperature Controller cum Indicator - 1°C resolution.
- ❖ Energy efficient tubular heater.



Optional Accessories:

- ❖ Gabled lid of stainless steel, specially designed so that the condensate drips on sides.
- ❖ A pack of float balls for containing evaporation loss. These balls will form a movable and flexible layer that facilitate easy removal and placement of test tubes, racks without the hassles of opening and closing a lid.
- ❖ Suitable test tube racks.

Model	Size L/W/D in cm	Racks (optional)		Optional SS Gabled Lid Model	Optional Float Ball Packs Model	Heater	App. Shipping Details Net / Gr. / Volume
		Kahn 30 #6205	Kahn 48 #6209				
#6701DI	27 x 13 x 12	-	1	#6901	#6990-50 x1	500 W	3.5 kg / 4.0 kg / 0.030 m ³
#6721DI	30 x 27 x 12	1	3	#6921	#6990-50 x 2	750 W	5.5 kg / 6.1 kg / 0.050 m ³

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request. Due to continuous improvement in design, the product supplied may have modified features.

Our Other Products: Autoclave Vertical: Standard / PAD / Fully Automatic / Jacketed, Autoclave Portable: Standard / Automatic, Autoclave SLL Series: Vacuum, Autoclave SLE series • Media Preparator • Table Top Autoclaves: Normal / Class B • UV Sterilization Cabinet • Rotary Vacuum Evaporator • Vacuum Controller • Chilled Water Circulator • PTFE Diaphragm Vacuum Pump • BOD Incubator • BSC Clean Bench • Water Still / Stone-fin / Automatic Water Still • Water Bath: Stirred (Constant) / Serological / Concentric Rings / Shaking / Refrigerated Circulating Bath • Round Bath: Clear Bottom / Immersion Series / ECO Series • Ultrasonic Bath: Digital / LCD • Dry Bath • Colony Counter • Hot Plate with Magnetic Stirrer • Anaerobic Jar: SS / Polycarbonate • Pipette & Petri Dish Sterilizing Box • Test Tube Rack • ESR Westergren & Wintrobe Apparatus • Slide Staining Stand / Rack / Trough • Slide Trays • Slide Cabinet



MEDICA INSTRUMENT MFG. CO.

105, Vasan Udyog Bhavan, Sun Mill Compound, Lower Parel (West)
Mumbai 400013, INDIA • ☎ +91-22/6618-9999 • 📠 +91-22/6618-9996
✉ sales@medicainstrument.com • www.medicainstrument.com



0015_1005-3-1