

'HANNA' MAGNETIC MINI-STIRRERS

Hanna HI 180 series is compact, lightweight, inexpensive, these stirrers are ideal for laboratories, schools and anywhere where space is a concern. Easy to handle, chemical resistant housing.

MODEL	CAPACITY	SPEED	BODY
HI 180	1ltr.	100Rpm to 1000Rpm	ABS round
HI 190M	1ltr.	100Rpm to 1000Rpm	ABS
HI 200M	1ltr.	100Rpm to 1000Rpm	S.S.cover

Complete with stir paddle and manual.



'HANNA' MAGNETIC STIRRERS HEAVY DUTY

Large capacity, heavy duty stirrers, can be used in lab as well as in production, with stainless steel cover, 220V AC mains operated. Complete with stir paddle and manual.

MODEL	CAPACITY	SPEED
HI 300N	2.5 ltr.	100 - 1000Rpm
HI 310N	5.0 ltr.	100 - 1000Rpm
HI 301N	2.5 ltr.	100-500 /100-1000Rpm
HI 311N	5.0 ltr.	100-500 /100-1000Rpm



'REMI' MAGNETIC STIRRERS

Stainless steel top housing series of magnetic stirrers are compact, with PMDC motors, lightweight, accurate speed control upto 1200Rpm, Ideal for laboratories, schools etc., Easy to handle.

MODEL	CAPACITY
MS-500	500ml.
ML 1	1 ltr.
ML 2	2 ltr.
ML 5	5 ltr.

Complete with stir paddle and manual.



'REMI' MAGNETIC STIRRERS WITH HOTPLATE

These are similar to ML series but with additional Stainless steel top hotplate. They are compact, with PMDC motors, lightweight, accurate speed control upto 1200Rpm, Ideal for laboratories, schools etc., Easy to handle. Specially designed for corrosive environment.

MODEL	CAPACITY	HEATING CAPACITY
MLH 1	1 ltr.	150W
MLH 2	2 ltr.	300W
MLH 5	5 ltr.	500W

Complete with stir paddle and manual.



'REMI' ROTAMANTLES

Specially designed to meet laboratory requirement of convenient stirring in flask with uniform heating by heating mantle. Accurate step less speed control up to 1200RPM,

MODEL	CAPACITY	HEATING CAPACITY
RML 1	1 ltr.	300W
RML 2	2 ltr.	400W
RML 3	5 ltr.	500W

