



VACUUM PUMP

'TIC' VACUUM PUMP CUM COMPRESSOR OIL FREE

Need of every lab for filtration under vacuum / pressure, for leak test, air/gas sampling, vacuum application etc.

- * Mini Oil free dry pump
- * Maintenance free
- * Built-in micro filter & vacuum and pressure regulators.

MODEL : TV 40(S)

Vacuum	: 23" of HG;
Pressure	: 35PSI(2BAR)
Free Air Delivery	: 35LPM;
Motor	: 1/3 HP
Watts	: 250W SINGLE PHASE,
Power	: 230V AC, 50Hz
Weight	: 7.5 Kg (Approx.)
Dimension	: 190mm x 150mm x 240mm (Approx.)



'TIC' VACUUM PUMP DOUBLE STAGE MODEL: TV-80

Vacuum	: 27" of Hg.
Free Air Delivery	: 45 LMP;
Motor	: 1/3 HP
Watts	: 250 W;
Phase	: Single 50 Hz
Power	: 230V AC ;
Weight	: 11 Kg (Approx.)
Dimension	: 190mm x 150mm x 340mm.



'TIC' VACUUM PUMP DOUBLE STAGE MODEL: TV-4025

Vacuum	: 27" of Hg.
Pressure	: 35 PSI; (2 BAR)
Motor	: 1/3 HP
Watts	: 250 W;
Phase	: Single 50 Hz
Power	: 230V AC
Weight	: 11 Kg (Approx.)
Dimension	: 190mm x 150mm x 340mm.

Complete with 2 gauges.



'TIC' VACUUM PUMP DOUBLE STAGE-DIRECT DRIVE

- * Rotary vane type pump. Alloy casting body.
- * Silent smooth working.
- * Built-in Air ballast, fitted with Crompton/Godrej Motor,
- * Operation on 220/230V AC, Mains 50 cycles. Single phase
- * Ultimate no load vacuum 0.001mm.Hg.

Supplied complete with 3 core wire, Vacuum Gauge, Moisture Trap and 1st filling vacuum oil.

Capacity / Min. : 100 Lit.

Capacity / Min. : 200 Lit.

Capacity / Min. : 300 Lit.

